

LOVIBOND® WATER ANALYSIS

TINTOMETER® GROUP



Lovibond® ShuddhCheck™

Your gateway to better Water Testing

INSTRUCTION GUIDE- IRON LR COLOUR CARD

WORKING RANGE: 0 – 1.5 mg/l (ppm)

MEASURING RANGE: 0, 0.2, 0.3, 0.4, 0.5, 1.0, 1.5 mg/l as Fe



Test Procedure

1. Take 1 ml water sample in a plastic vial and mix it with 2ml of DI water
2. Add 12 drops of Iron LR Reagent 1 and mix well
3. Wait for 2 minutes
4. Add 1 drop of Iron LR Reagent 2 and mix well
5. Wait for 2 minutes
6. Notice the colour change of the liquid surface and compare it with the colour card provided with the kit
7. The matching colour corresponds to the reading of Iron in mg/l (ppm) in the water sample

This will give the Iron value directly in mg/l (ppm) as Fe in that water sample.

No. of Tests: 100

Order Code: IN415070MK

Reorder Code: IN415071MKR

Tintometer India Private Limited

Door No. 7-2-C-14, 2nd, 3rd & 4th Floor

Sanathnagar Industrial Estate | Hyderabad-500018

Telangana | India



Notes

- * Cleaning: Thoroughly rinse out sample container after use.
- ** Test reagents are for chemical testing only. Not meant for human consumption; Keep out of reach of children. For industrial use only.



For more information please reach out to us:

☎ 1800 599 3891/92 | ✉ indiaoffice@lovibond.in | 🌐 www.lovibond.in

🌐 <https://in.linkedin.com/company/tintometer-india-private-limited>