

LOVIBOND® WATER ANALYSIS

TINTOMETER® GROUP



## Lovibond® ShuddhCheck™

Your gateway to better Water Testing

### INSTRUCTION GUIDE- IRON LR- I COLOUR CARD

WORKING RANGE: 0 – 2 mg/l (ppm)

MEASURING RANGE: 0, 0.2, 0.3, 0.6, 0.8, 1.0, 1.5, 2.0 mg/l as Fe



### Test Procedure

1. Take 10 ml water sample in a glass vial
2. Add 12 drops of Iron LR Reagent 1 and mix well
3. Add 2 drops of Iron LR Reagent 3 and mix well
4. Wait for 2 minutes
5. Notice the colour change and compare it with the colour card provided with the kit
6. The matching colour corresponds to the reading of Iron in mg/l (ppm) in the water sample  
This will give the Iron value directly in mg/l (ppm) as Fe in that water sample.

No. of Tests: 100

Order Code: IN415310MK

Reorder Code: IN415311MKR

## Tintometer India Private Limited

Door No. 7-2-C-14, 2nd, 3rd & 4th Floor

Sanathnagar Industrial Estate | Hyderabad-500018

Telangana | India



## Notes

- \* Cleaning: Thoroughly rinse out sample container after use.
- \*\* Test reagents are for chemical testing only. Not meant for human consumption; Keep out of reach of children. For industrial use only.



**For more information please reach out to us:**

☎ 1800 599 3891/92 | ✉ [indiaoffice@lovibond.in](mailto:indiaoffice@lovibond.in) | 🌐 [www.lovibond.in](http://www.lovibond.in)

🌐 <https://in.linkedin.com/company/tintometer-india-private-limited>